


Dental matrices and wedges



ILO's:

By the end of the lecture the student should be able to:

Describe about Objectives and ideal requirements of matricing

Explain the various forms of matricing

Describe the purposes and type of wedges


Explain about various wedging techniques

Class- II





Class- III or IV

Not using Matrix & Wedge




Using Matrix & Wedge



Why we need Proper Contacts and Contours

- Prevents the food from impinging between the teeth
- Prevent periodontal problems as well as formation of dental caries
- Stabilizes the dental arches
- Distributes the forces of mastication.



What is matricing?

Matricing is a procedure wherein a temporary wall is created opposite to axial walls, surrounding areas of tooth structure that were lost during preparation



The primary objective of the matrix is to restore anatomic contours and contact areas.

The Ideal requirements of a matrix include :

(1) Rigidity

(2) **Versatility** (The matrix band should be able to adapt to almost any size and shape of tooth)




(3) Establishment of proper anatomic contour

(4) Restoration of correct proximal contact

(5) Prevention of gingival excess(Overhanging)

(6) Ease of application

(7) Ease of removal




Sturdevant's Art and Science of Operative Dentistry.
6th Edition. Chapter 14, Page 394-403
5th Edition. Chapter 17, Page 760-773

MATRIX SYSTEMS

1. Universal Matrix (Tofflemire Matrix)



2. Pre-contoured matrix (Sectional matrix)



3. Auto-matrix



4. Celluloid / Polyester strips



UNIVERSAL MATRIX (TOFFLEMIRE)

BY, B.R.TOFFLEMIRE

✓ Commonly used to restore Class II preparations.

2 parts

Matrix Retainer - It is an instrument used to hold matrix band in position


Matrix band

✓ **Disadvantage** - Flat band must be contoured.



CE

Universal Matrix Retainer Tofflemire



PREMIUM
INSTRUMENTS



<https://www.youtube.com/watch?v=upzeWfgvpQ8>

UNIVERSAL MATRIX



CONTOURING THE BAND (To achieve proper convex contour and proper adaptation)



DO IT ON PAPER PAD



USE EGG BURNISHER



POSITIONING OF THE UNIVERSAL MATRIX BAND IN THE RETAINER



PLACE THE GROOVE GINGIVALLY




LOOSEN THE SMALL KNARLED NUT



POSITIONING OF THE UNIVERSAL MATRIX BAND IN THE RETAINER




The convex surface (upward) & one concave (downward). The convexity of matrix band should put on the occlusal surface.
Form a loop




**PASS IT IN THE GINGIVAL SLOT AND
TIGHTEN THE LARGE KNARLED NUT.**


POSITIONING OF THE UNIVERSAL MATRIX BAND IN THE RETAINER



**FORM LOOP WITH BAND AND
POSITION**



POSTIONING THE MATRIX RETAINER AND THE BAND OVER THE TOOTH




The matrix band should be 1 mm higher than the adjacent tooth marginal ridge and 1 mm below the gingival contact area

DO NOT TRAP THE RUBBER DAM IN THE BAND





No gap between the band and the tooth

Use explorer tip (with pressure) between matrix band and gingival seat to ensure proper adaptation of the band to gingival margin.



IN CASE OF DEEP PROXIMAL MARGINS(SUBGINGIVAL RESTORATIONS) AND SHALLOW MARGINS



USE THE MOUTH MIRROR FACIALLY AND LINGUALLY TO EVALUATE THE PROXIMAL CONTOUR



PRE CONTOURED MATRIX (SECTIONAL MATRIX)

PALODENT

1. The Sectional matrix is the best way to achieve strong contact point in a class II composite restorations.
2. Glass-fiber reinforced plastic tines maintain a fit on the tooth that complements the wedges and minimizes flash.
3. Ring design promotes exceptional stability in the forceps during placement. The function of the ring is to adapt the sectional matrix in the cervical area.
4. V-shaped tines allow you to wedge from both sides before or after ring placement.



Sectional Matrix (PALODENT/TRIDENT systems)



AUTOMATRIX



AUTOMATRIX



ADVANTAGE –

1. Autolock loop can be positioned on the facial and lingual surface with equal ease.
2. Improved visibility , because of absence of retainer

DIASDVANTAGE –

1. The bands are not precontoured
2. Development of physiological proximal contours is difficult.



POLYESTER / CELLULOID MATRIX SYSTEMS


Used in class III and class IV preparations



MYLAR MATRIX STRIPS



To restore Class-III and IV Cavity Preparation




WEDGES

Wedges are the devices that are placed to create rapid separation during tooth preparation and restoration and are placed in the gingival embrasure area

FUNCTIONS OF A WEDGE

1. **Rapid separation of the tooth.**
2. **Protecting the underlying gingiva from unexpected trauma**
3. **To prevent gingival overhangs of restoration.**
4. **Compensating for the thickness of the matrix band to achieve proximal contact.**





For posterior restorations, the wedge is placed from the larger to the smaller interproximal gingival embrasure, typically from a lingual approach, just apical to the gingival margin.



TYPES OF WEDGES

Wooden wedges



Plastic wedges

FENDERWEDGE®
Interproximal Tooth Shields



Garrison

CLASSIFICATION OF WEDGES

Based on the method of fabrication

Custom made

Prefabricated

Based on the material used

Wooden wedges

Resin /Plastic wedges

WODDEN WEDGES



Advantage





Number designates wedge length in mm (13=13mm)

- Can be easily customized (Can be cut and trimmed by scalpel blade)
- They absorb water interiorly and swells up, which presses against the matrix to give more retention

LIGHT TRANSMITTING WEDGES FOR COMPOSITE RESTORATIONS –

Transparent plastic wedges which are available with light reflecting property, helps to assist in directing the light into the interproximal areas during curing



PLASTIC WEDGES



Plastic wedges allows the visible light to pass through helpful in Composite restorations

PLASTIC WEDGES




The Wand applicator provides the simplest placement technique ever.

No instrument such as cotton plier is required to place the wedge

Wedges with handles

WEDGE WANDS


PLASTIC WEDGES



G-Wedge™



Perfect wedge, perfect contact.



Traditional wedge placement with
tweezer or cotton plier

G-Wedge™ gives you the same basic
wedge shape as the **Wedge**
Wands® without the applicator wand.

DIAMOND ANATOMICAL WEDGE



Diamond anatomical Wedge requires less effort to place than traditional wedges due to the diamond shaped cut-out in the tip, which collapses during insertion. The flexibility of the Diamond Wedge creates an excellent marginal adaptation.

FENDER MATRIX AND WEDGE



FenderWedge is a combination of a steel band and a plastic wedge.

Fender Wedge® will lend you a hand by "pre-separating" the teeth while you prepare the tooth.

Protect the adjacent teeth

Large: 2.3mm

Medium: 1.9mm

Small: 1.4mm

X-Small: 1.0mm





As a added bonus, the firm stainless steel shield protects the adjacent tooth from inadvertent contact with the bur.

PLASTIC WEDGES

COMPOSI –TIGHT 3D FUSION ULTRA ADAPTIVE WEDGE




A soft, rubberized material is fused to the firm plastic core.

Advantages

It prevents backout of the wedge


The wedge gently molds itself to irregularities as you insert it providing an unsurpassed seal.

WEDGE PLACEMENT



If the wedge is placed occlusal to the gingival seat, the band will be pressed into the preparation, creating an abnormal concavity in the proximal surface of the restoration.

It should not interfere with the curvature of the matrix band




The wedge should be apical to the gingival seat that the band will be held tightly against the gingival margin.

MODIFIED WEDGING TECHNIQUES



I. Piggy back wedging

In this technique, smaller wedge is placed over the first wedge, This method is useful when proximal box is shallow gingivally or interproximal gingival tissue level recession





2. Double Wedging

Double wedging, as the name indicates one wedge from facial and other in lingual embrasure, may be used with buccolingually wide proximal boxes to provide maximal closure of band along gingival margin.




3. Wedge Wedging

Another technique may be used on the mesial aspect of maxillary first premolars to adapt the matrix to the fluted (i.e., Concave) area of the gingival margin ,*second wedge inserted between the first wedge and matrix band* .



Fluting results in opening between matrix and gingival margin



Wedge-wedging
a matrix, ready for compound;
wedges inserted from lingual or
facial embrasure, whichever is larger

REMOVAL OF THE MATRIX AND WEDGE



First remove the retainer, then
matrix band followed by wedge

DISPOSAL OF MATRIX BAND OR ANY SHARP ITEMS

Needles, scalpel blades, Local anesthetic cartridges, Burs, Matrix bands must be disposed of in the yellow sharps container



Reference:

1. Sturdevant's Art and Science of Operative Dentistry. 6th Edition
Chapter 14, Page 394-403
2. Sturdevant's Art and Science of Operative Dentistry. 5th Edition
Chapter 17, Page 760-773
3. <https://www.garrisondental.com/en/wedges>
4. <https://www.dentsplysirona.com/en-ca/products/restorative/accessories/matrix-systems.html>